




International Symposium for Thin Layer Chromatography

Helsinki, 11th – 13th June 2008

**HPTLC/ DESI- MS for the analysis of proteins and
peptides**

Michael Schulz
PLS R&D LSS SIS

HPTLC/ DESI- MS of peptides from tryptic digest



HPTLC/ DESI- MS of peptides from tryptic digest



Tryptic digest, HPTLC and staining performed by Merck KGaA

HPTLC/ DESI- MS and HPLC/ ESI- MS performed by Organic and Biological Mass Spectrometry Group, Chemical Science Division, Oak Ridge National Laboratory, Oak Ridge

Lit.: Using HPTLC/ DESI- MS for peptide identification in 1D separations of tryptic protein digests, Anal Bioanal Chem (2008) 391: 317-324

Experimental / HPTLC



Stationary phases:

ProteoChrom® HPTLC Silica gel F254s, from Merck KGaA

Format 20x10 cm

Layer thickness 100 µm

Special binder composition

ProteoChrom® HPTLC Cellulose, from Merck KGaA

Format 10x10 cm

Layer thickness 100 µm

High layer density

Mobile phases:

2-butanol/pyridine/acetic acid/ water (30/20/6/24) (v/v/v/v) and

2-butanol/pyridine/ammonia/ water (39/34/10/26) (v/v/v/v)

Instruments:

ATS4, Normal flat bottom chamber, TLC plate heater and DigiStore 2 from CAMAG (Muttens, Switzerland)

DS 20 from DESAGA (Nürnberg, Germany)

Experimental / HPTLC



Staining:



Ninhydrin: The plate is sprayed with ninhydrin solution and heated up to 120°C for 2 min



ProteoChrom® Color: The sheet is sprayed with ProteoChrom® Color Reagent. After drying at room temperature the sheet is sprayed with ninhydrin solution and heated up to 120°C for 1-2 min

Experimental / HPLC



Instrument: Agilent 1100 Series Capillary LC System
Column: Acclaim PepMap 100 C18 column (1 mm x 150 mm; 3 µm particle size; Applied Biosystems)
Mobile phase: water (0.1% formic acid) (A) and 80% acetonitrile/20% water (0.1% formic acid) (B). (Gradient elution: 0 min with 96% A : 4% B; 30 min to 80% A : 20 % B; 20 min 80% A : 20 % B).

Experimental / MS



Instruments

MS: LCQ Deca ion trap or

LTQFT Ultra Hybrid Mass Spectrometer (Thermo Scientific, San Jose, CA, USA)

DESI Spray emitter: MicrolonSpray II (MDS Sciex, Concord, Ontario, Canada)

DESI conditions

Solvent: H₂O

Solvent flow rate: 5 µL/min


Spray voltage: 4 kV

Gas flow (N₂): 1,7 l/h


Scan rate: 100 µm/s

AGC on

Experimental / MS



Experimental / MS




1D- HPTLC/ DESI- MS



ProteoChrom® Cellulose

Lys Myo Cyt Cas BSA




2-butanol/pyridine/acetic acid/water
(30/20/6/24) (v/v/v/v)

ProteoChrom® Color

ProteoChrom® Silica gel 60 F254s

Lys Myo Cyt Cas BSA



2-butanol/pyridine/ammonia/water
(39/34/10/26) (v/v/v/v)


Ninhydrin

1D- HPTLC/ DESI- MS




Tryptic digest of Cytochrome C separated on:

ProteoChrom® Cellulose



ProteoChrom® Silica gel 60 F254s



Mass spectrometer: LCQ Deca ion trap
Application volume: 7 μ l
Concentration: 2 mg/ml
Total amount of protein: ~ 14 μ g

1D- HPTLC/ DESI- MS



Tryptic digest of Cytochrome C separated on:

ProteoChrom® Cellulose

GDVEK**GK**KIF VQKCAQCHTV
EKGGK**HK**TGP NLHGLFGRKT
GQAPGFSYTD ANKNK**GITWG**
EETLMEYLEN **PKK**YIPGTKM
IFAGIK**K**KGE REDLIAYLKK ATNE

Sequence coverage 72%

ProteoChrom® Silica gel 60 F254s

GDVEK**GK**IF VQKCAQCHTV
EKGGK**HK**TGPNLHGLFGRKT
GQAPGFSYTD ANKNK**GITWG**
EETLMEYLEN **PKK**YIPGTKM
IFAGIK**K**KGE REDLIAYLKK ATNE


Sequence coverage 58%

1D- HPTLC/ DESI- MS




Tryptic digest of Myoglobin separated on:

ProteoChrom® Cellulose



ProteoChrom® Silica gel 60 F254s



Mass spectrometer:	LCQ Deca ion trap
Application volume:	7 μ l
Concentration:	2 mg/ml
Total amount of protein:	~ 14 μ g

1D- HPTLC/ DESI- MS



Tryptic digest of Myoglobin separated on:

ProteoChrom® Cellulose

GLSDGEWQQV LNVWGKVEAD IAGHGQEVLI
RLFTGHPETL EKFDKFKHLK TEAEMKASED
KKHGTVVLT ALGGILKKKG HHEAELKPLA
QSHATKHKIP IKYLEFISDA IHVLHSHKHP
GDFGADAQGA MTKALELFRN DIAAKYKELG
FQG

Sequence coverage 68%

ProteoChrom® Silica gel 60 F254s

GLSDGEWQQV LNVWGKVEAD IAGHGQEVLI
RLFTGHPETL EKFDKFKHLK TEAEMKASED
LKKHGTVVLT ALGGILKKKG HHEAELKPLA
QSHATKHKIP IKYLEFISDA IHVLHSHKHP
GDFGADAQGA MTKALELFRN DIAAKYKELG
FQG


Sequence coverage 62%

Sequence coverage 86%


1D- HPTLC/ DESI- MS



ProteoChrom® Cellulose




ProteoChrom® Silica gel 60 F254s



- Separation of peptides on Silica gel and Cellulose mainly according hydrophilicity

2D- HPTLC/ DESI- MS




2D- HPTLC/ DESI- MS





Tryptic digest of Cytochrome C separated on:

ProteoChrom® Cellulose



Application volume: 10 μ l
Concentration: 2 mg/ml
Total amount of protein: ~ 20 μ g

Chemical image



LCQ Deca ion
trap Mass
Spectrometer

LTQFT Ultra
Hybrid Mass
Spectrometer

2D- HPTLC/ DESI- MS



		HPLC/ ESI- MS	1D-HPTLC/ DESI-MS	2D-HPTLC/ DESI-MS
ProteoChrom® HPTLC Silica gel 60 F254s	Cytochrome C		58%	
ProteoChrom® Cellulose	Cytochrome C		72%	81%
Acclaim PepMap 100 C18	Cytochrome C	92%		
ProteoChrom® HPTLC Silica gel 60 F254s	Myoglobin		62%	
ProteoChrom® Cellulose	Myoglobin		68%	74%
Acclaim PepMap 100 C18	Myoglobin	84%		

Conclusions



- Limited separation efficiency of 2D- HPTLC
- Low DESI signal for peptides with low hRf- value
- + Direct sampling from the plate under ambient conditions
- + Lane scan and spot sampling mode possible
- + 2D-HPTLC/ DESI- MS Sequence Coverage is comparable to HPLC/ ESI- MS
- + Separation can be stored on the HPTLC plates
- + Numerous separations with every HPTLC run

Appreciation



- Merck KGaA
Susanne Schorcht
- Organic and Biological Mass Spectrometry
Group, Oak Ridge National Laboratory
Sofie P. Pasilis
Vilmos Kertesz
Gary J. Van Berkel