


HPTLC

For the analysis of botanical materials and medicinal plants

Eike Reich

CAMAG Laboratory
Sonnenmattstrasse 11
4132 Muttenz / Switzerland



What will be discussed?

1. What exactly is HPTLC?
2. The elements of a standardized HPTLC methodology
3. Changing TLC methods into HPTLC methods
4. Development and validation of ID Methods
5. How new developments in HPTLC may affect the analysis of plants

Introduction

- ▶ Thin-layer chromatography (TLC) always was and still remains an important tool for the analysis of plants.
- ▶ Today there are two principal applications in this context: research and quality control. Both benefit from the advantages of the planar off-line principle and also in particular from low cost, simplicity, and flexibility.
- ▶ For decades TLC is integral part of monographs for medicinal plants in all pharmacopoeias and the primary method of identification.
- ▶ Growing expectation regarding performance characteristics have brought TLC methods to the limits.
- ▶ Since the turn of the century pharmacopoeias recognize the technical progress in instrumentation and improvements offered by high performance plates.
- ▶ Most recently HPTLC is being discussed as alternative to classical TLC.

The planar off-line principle

application

development

detection and quantitative evaluation

documentation



FLEXIBILITY

But what about reproducibility?



The central problem:

- ▶ When two (or more) labs do the „same“ or think that they are doing the same, the results are not necessarily equal.
- ▶ Reason is the general method description in the Pharmacopoeias (EP 2.2.27, USP <201>, <621>, PhPRC ap. VI) which define „suitable“ equipment and give ranges instead of values.
- ▶ JPXV 2.03 is still centered around self made 20x20 TLC plates !
- ▶ A table (EP) or a result description (USP, JP) can only define the most important aspects of a TLC chromatogram. That leaves room for interpretation. An atlas (PhPRC) provides clear guidance.
- ▶ Example: how can a color be described correctly?

Most text books are even worse ...

6


What is blue?




PhEur: Possible choices in methodology

- ▶ TLC layer
- ▶ Manual application
- ▶ Transparent container (Pickle jar?)
- ▶ UV-Lampe (λ ?)
- ▶ Manual spraying / immersion
- ▶ HPTLC layer
- ▶ Automatic application
- ▶ Automatic Developing Chamber
- ▶ Scanner
- ▶ Automatic immersion / spraying



Identification of *Acanthopanax*



Identification of Peonies



Identification of Fleece flower



11





Goal:

International standardization of HPTLC


- ▶ What
 - HPTLC definition
 - Methodology
 - Equipment
- ▶ Why
 - Reproducibility of results
 - Validity of official methods
 - Quality assurance in a globalized world
 - Quality of published research
- ▶ How
 - International collaboration
 - Top down
 - Publication

TLC or HPTLC

- ▶ Pharmacopoeias see difference primarily in the plate yet assume similar results




Comparison TLC-HPTLC of flavonoids





← TLC plate 20 x 20 cm
(135 mm)

HPTLC plate 20 x 10 cm
(60 mm) ↓



14

TLC or HPTLC?



What is TLC?

- ▶ Chromatography for the poor (cheap)
- ▶ Simple manual chromatography for everyone (students?)
- ▶ Rapid
- ▶ Flexible
- ▶ Reference and test solution side by side
- ▶ “Just” qualitative, preliminary estimation at best

- ▶ Unpredictable
- ▶ Unreliable

→ Manual technique, simple instruments, TLC plates

What is HPTLC?

High Performance Thin-Layer Chromatography

TLC for the 21st century

- ▶ Instrumental TLC
 - Application
 - Development
 - Documentation
 - Densitometry
- ▶ Truly „plug and play“
- ▶ Fully cGMP compliant

A new concept

- ▶ Suitable instruments
- ▶ Scientific basis
- ▶ Standardized methodology
- ▶ Validated methods

“HPTLC” on TLC plates?

- ▶ What is the point?
 - Saturation 1h vs. 20 min
 - Twice (10x) the solvent volume
 - 3x the developing time (15 vs. 6 cm)
 - Same cost per plate (20x20 vs. 10x10 cm)

yet

 - Less resolution
 - No control of the development process

Standardization of methodology

- ▶ Plate setup and handling
- ▶ Sample application (as band)
- ▶ Chamber geometry and saturation
- ▶ Humidity control
- ▶ Developing distance
- ▶ Derivatization procedure
- ▶ Documentation (electronic images)
- ▶ Evaluation



SOP

A Standardized Approach to Modern High-Performance Thin-Layer Chromatography (HPTLC). Reich, E., Schibli, A., (2004) J. Planar Chromatogr. 17, 438-443

SOP for HPTLC

- ▶ Should be the basis for all work (in participating labs)
- ▶ Applies to all methods
- ▶ All deviations need to be recorded
- ▶ Our SOP is in full compliance with PhEur, USP, ChP

Available at: www.camag-laboratory.com (homepage)


The basic HPTLC setup

Application

Development

Derivatization

Documentation




Standardized basic equipment

- ▶ Clear specification of HPTLC plates (type, format, manufacturer)
- ▶ Software control for “absolute” reproducibility of all parameters
- ▶ Independence from environmental factors (humidity !)
- ▶ Eliminating human factors

- ▶ Emulation by “manual” operation still possible to a certain degree

Sources of methods

- ▶ European Pharmacopoeia (EP)
 - New monographs feature TLC and HPTLC in parallel
- ▶ British, French, German, Swiss Pharmacopoeias
 - Offer specific monographs not found in EP
- ▶ The USP Dietary Supplement Compendium
 - TLC and state of the art HPTLC
- ▶ Chinese Pharmacopoeia
 - HPTLC atlas as a supplement of 2005 ed.
- ▶ Japanese Pharmacopoeia
 - Only simple TLC



More sources of methods

- ▶ American Herbal Pharmacopoeia (19 monographs)
- ▶ Quality Standards of Indian Medicinal Plants
(8 volumes)
- ▶ Indian Herbal Pharmacopoeia (54 monographs)
- ▶ Quality Standards of Traditional Chinese Medicines
(Chinese only)
- ▶ Wagner, H. and Bladt, S. „Plant Drug Analysis“

Converting existing TLC methods to HPTLC

Assumptions:

- ▶ Results on HPTLC and TLC plates are similar if
 - No changes are made to chromatographic system
 - Same equipment is used
 - Original TLC method was optimized
- ▶ Due to the higher separation power HPTLC plates
 - Usually give improved result
 - Require shorter developing distance → less time
- ▶ HPTLC with no instruments looks bad, but so does TLC

Converting existing TLC methods to HPTLC

Practical aspects (I)

- ▶ Do not change chromatographic system (chamber configuration, mobile phase, stationary phase)
- ▶ Reduce application volume (generally) to 1/5 (typically 2 μ L)
- ▶ Employ standardized methodology (based on SOP):
 - Fixed (x, y) application positions, (e.g.) 8 mm bands
 - Use 60 mm developing distance
 - Fixed drying time and temperature
 - Use dipping instead of spraying if possible
 - Fixed waiting times between derivatization and evaluation
 - Obtain multiple images (if possible) of plate (e.g. UV 254, UV 366 prior to derivatization, and white light, UV 366 after derivatization)

Converting existing TLC methods to HPTLC


Practical aspects (II)

- ▶ Use same samples on TLC and HPTLC
- ▶ Evaluate whether changes in the result can
 - Fall under “additional weak zones may be seen”
 - Are due to natural variability of plant
- ▶ Color description is always subjective: e.g.
 - blue, bluish
 - blue white
 - blue green, etc.
- ▶ R_f is predictable only for validated methods!

Example: USP method for chamomile



TLC



HPTLC
application volume 2 μ L

Developing an ID method from scratch

- ▶ Review literature for related plants
- ▶ Obtain multiple samples form different accessions
- ▶ Obtain samples of related plants and known adulterants
- ▶ Optimize sample preparation and detection
- ▶ Avoid toxic solvents
- ▶ Start with silica gel, select mobile phase

Validation of qualitative methods




Validation of high-performance thin-layer chromatographic methods for the identification of botanicals in a cGMP environment. Reich, E., Schibli, A., DeBatt, A. (2008) J. AOAC Int. 91, 13-20.

Envisioning the Future...


Methods for identification of plants

- ▶ SOP is the basis for an HPTLC method template (instruments)
 - [software](#) template
- ▶ A (PhEur, USP, etc.) monograph or validated methods is basis for HPTLC method document → Method for identification of [Thyme](#)
- ▶ Result table/description is replaced by a reference image
- ▶ Corresponding instrument methods are derived from software template
 - [instrument](#) method
- ▶ Results on plate are qualified by a System Suitability Test (SST)
- ▶ Comparison of results with reference images

Thyme leaf plates 1 - 3



Thyme leaf



34



Goal:

The CAMAG HPTLC method collection

- ▶ Identification of 150 (+) plants from East and West
- ▶ Predictable results everywhere
- ▶ Convenient comparison against reference images
- ▶ Compatible with (current and future) harmonized description of HPTLC
- ▶ Based on ATS4, ADC2, Visualizer and Software

35



Goal:

CAMAG Laboratory Network

- ▶ Collaboration on development and validation of HPTLC methods for identification of plants
- ▶ Setting global quality standards for HPTLC
- ▶ Providing training and research opportunities
- ▶ Offering analytical services for customers

Who is part of the network?


- ▶ CAMAG Laboratory Muttenz
- ▶ CAMAG Scientific Inc. Wilmington (NC)
- ▶ Shanghai University of TCM (Prof. Wang)
- ▶ University of Applied Sciences Wädenswil (Prof. Meier)
- ▶ University of Regensburg (Prof. Heilmann)
- ▶ University of Barcelona (Prof. Canigüeral)
- ▶ University Sapienza Rome (Prof. Nicoletti)
- ▶ University of Graz (Prof. Bauer)

Examples from other labs (HS Wädenswil)

Oswego tea

4. Results

Fig.1) NP reagent and PEG reagent, UV 366 nm



Track	Volume	Sample	Track	Volume	Sample
1	2 µL	Oswego tea; herb 1	9	2 / 5 µL	Chlorogenic acid, rosmarinic acid (with increasing R _f)
2	4 µL	Oswego tea; herb 1	10	4 µL	Oswego tea; herb 1
3	6 µL	Oswego tea; herb 1	11	4 µL	Oswego tea; herb 2
4	2 µL	Oswego tea; herb with flower	12	4 µL	Oswego tea; herb 3
5	4 µL	Oswego tea; herb with flower	13	4 µL	Oswego tea; herb 4
6	6 µL	Oswego tea; herb with flower	14	4 µL	Oswego tea; herb with flower
7	2 µL	Rutin, narirutin, hyperoside, quercetin (with increasing R _f)	15	4 µL	Oswego tea; flower
8	2 µL	Naringin, naringenin (with increasing R _f)			

System suitability test

Rutin: orange fluorescent zone at R_f ~ 0.32 (UV 366 nm).


Hyperoside: orange fluorescent zone at R_f ~ 0.50 (UV 366 nm).

Examples from other labs (Uni Regensburg)

Coptis rhiz.

4. Results

Fig. 1) UV 366nm



Track	Volume	Sample	Track	Volume	Sample
1	10 µL	Palmatine and berberine	9	10 µL	Coptis teeta rhizome
2	10 µL	Coptisine	10	10 µL	Chinese corktree bark
3	10 µL	Coptis rhizome # 1	11	10 µL	Chinese mahonia bark
4	10 µL	Coptis rhizome # 2	12	10 µL	Tinospora root
5	10 µL	Coptis rhizome # 3	13		blank
6	10 µL	Coptis rhizome # 4	14		blank
7	10 µL	Coptis rhizome # 5	15		blank
8	10 µL	Coptis deltoidea rhizome			

System suitability test

Palmatine: fluorescent zone at $R_f \sim 0.14$.

Berberine: fluorescent zone at $R_f \sim 0.23$.

Coptisine: fluorescent zone at $R_f \sim 0.60$.



Quick access and search

P340 P340 Probensammlung Test 7 Flavonoide



Search...

Methods (1)
Plate Analysis (28)

Today

- Today
- Yesterday
- A few days ago
- Earlier this month
- Earlier this year
- A long time ago

File Preview



Tr.	Vial ID	Description	Vol.
1	R6424	Rosmarinic acid	2.0
+	R7637	Rutin	2.0
+	R7638	Hyperosid	2.0
2	S2930	Serpylli herba con...	2.0
3	S6522	Thymus serpyllum	2.0
4	S8797	Thymus serpyllum	2.0
5	S8798	Thymus serpyllum...	2.0
6	S6587	Thymian	2.0
7	S8457	Thymus zygis Loef...	2.0
8	S8460	Thymus hyemalis...	2.0
9	S8461	Thymus moroderi...	2.0
10	S8462	Thymbra capitata (...)	2.0
11	S6097	Feldthymian	2.0
12	R6424	Rosmarinic acid	2.0
+	R7637	Rutin	2.0
+	R7638	Hyperosid	2.0

Name	Created by	Modified
Probensammlung P340	Rebekka Ambuehl	01.03.2011 12:04
P340_110302_02	Rebekka Ambuehl	09.05.2011 11:54
P340_110302_03	Rebekka Ambuehl	27.04.2011 08:40
P340_110302_04	Rebekka Ambuehl	15.04.2011 17:10
P340_110302_05	Rebekka Ambuehl	09.05.2011 11:54
P340_110302_06	Rebekka Ambuehl	16.06.2011 15:05
P340_110302_07	Rebekka Ambuehl	16.06.2011 15:09
p340_110302_08	Rebekka Ambuehl	08.04.2011 15:32
P340_110303_01	Rebekka Ambuehl	08.04.2011 15:36
P340_110303_02	Rebekka Ambuehl	08.04.2011 15:38
P340_110303_03	Rebekka Ambuehl	08.04.2011 15:40
P340_110303_04	Rebekka Ambuehl	08.04.2011 15:43
P340_110303_05	Rebekka Ambuehl	08.04.2011 15:45
P340_110303_06	Rebekka Ambuehl	08.04.2011 15:47
P340_110303_07	Rebekka Ambuehl	08.04.2011 15:50
P340_110303_08	Rebekka Ambuehl	08.04.2011 15:52
P340_110304_01	Rebekka Ambuehl	08.04.2011 15:54
P340_110304_02	Rebekka Ambuehl	15.04.2011 16:16
P340_110304_03	Rebekka Ambuehl	16.06.2011 15:05
P340_110304_04	Rebekka Ambuehl	15.04.2011 16:27
P340_110304_05	Rebekka Ambuehl	11.05.2011 10:34
P340_110304_06	Rebekka Ambuehl	09.05.2011 13:22
P340_110304_07	Rebekka Ambuehl	28.04.2011 11:43
P340_110304_08	Rebekka Ambuehl	15.04.2011 16:47
P340_110310_01	Rebekka Ambuehl	15.04.2011 17:17
P340_110310_02	Rebekka Ambuehl	11.05.2011 11:29
P340_110310_03	Rebekka Ambuehl	16.06.2011 15:11
P340_110310_04	Rebekka Ambuehl	09.05.2011 11:42
P340_110504_01_EC	Rebekka Ambuehl	11.05.2011 10:16
All fingerprints	Rebekka Ambuehl	29.06.2011 09:35
Clones	Rebekka Ambuehl	10.03.2011 09:06
Sageleaf	Rebekka Ambuehl	27.04.2011 09:48
Sageleaf	Rebekka Ambuehl	13.04.2011 16:36
schöneStandards	Rebekka Ambuehl	29.06.2011 09:36

A new software concept!

Explorer P340_110504_01_EC * x

Sequence Table Definition

Sequence table

Track	Vial ID	Description	Volume	Position	Type
1	R6424	Rosmarinic acid	2.0	A1	Standard
+	R7637	Rutin	2.0	A2	Standard
+	R7638	Hyperosid	2.0	A3	Standard
2	S2930	Serpylli herba concisa	2.0	B1	Sample
3	S6522	Thymus serpyllum	2.0	B2	Sample
4	S8797	Thymus serpyllum	2.0	B3	Sample
5	S8798	Thymus serpyllum pulvis	2.0	B4	Sample
6	S6587	Thymian	2.0	B5	Sample
7	S8457	Thymus zygis Loeffl. Ex L.	2.0	B6	Sample
8	S8460	Thymus hyemalis Lange	2.0	B7	Sample

Plate Layout Preview

Layout

Left

Center

Notes

OK Cancel

Plate Layout

Parameters

Plate Manufacturer **Merck**

Stationary phase **HPTLC plates silica gel 60 F**

ADC2 development Settings

Parameters

Mobile phase **Ethyl acetate, dichloromethane, formic acid, acetic ac**

Development type **Saturat**

Saturation (default)

Activation (default)

Drying (default)

Preconditioning

Notes

Visualizer Documentation Settings

Parameters

Image quality **Enhanced** HDRI **Restore default**

White R ▶

White RT ▶

White T ▶

254 nm ▶

366 nm ▶

Image sequence ▶

(Inactive)

Notes

OK

Cancel

Application

Development

Derivatization

Data Acquisition




ATS 4



ADC 2




Derivatization dip



Visualizer



Linomat 5



Chamber




Derivatization spray



Nanomat

TLC Steps



Finish Step Definition


Data Type



Overview

Chronological Illumination



Postderivatized 2a (2)



Postderivatized 1a (2)



Developed 1a (3)






Data Type



Overview

Tr.	Vial ID	Description
1	R6424	Rosmarinic acid
	+ R7637	Rutin
	+ R7638	Hyperosid
2	S2930	Serpylli herba concisa
3	S6522	Thymus serpyllum
4	S8797	Thymus serpyllum
5	S8798	Thymus serpyllum pulvis
6	S6587	Thymian
7	S8457	Thymus zygis Loeffl. Ex L.
8	S8460	Thymus hyemalis Lange
9	S8461	Thymus moroderi Pau ex M...
10	S8462	Thymbra capitata (L) Cav.
11	S6097	Feldthymian
12	R6424	Rosmarinic acid
	+ R7637	Rutin
	+ R7638	Hyperosid



Set

Batch set:

Selected Unselected


References

Mixture




Samples

S2930 Serpylli herba concisa
S8798 Thymus serpyllum pulvis
S8457 Thymus zygis Loeffl. Ex L.
S8461 Thymus moroderi Pau ex
S8462 Thymra capitata (L.) Ga
S6097 Feldthymian





Information

Three channels per track



PCA (Master thesis R. Ambühl, Uni Basel)



HPTLC as research tool

- ▶ Evaluation of column fractions
- ▶ Identification by MS
- ▶ Screening for bio-activity

Application of HPTLC-MS for the identification of flavonoids in herbal extracts

Bauer R¹, Meier M¹, Pferschy-Wenzig E¹, Woelkart K¹, Reich E²
¹ Institute of Pharmaceutical Sciences, Med-Faculty-University Graz, Universitätsplatz 4, 8010 Graz, Austria
² CAMMAG Laboratory, Sonnenmattstrasse 11, 4132 Mültenz, Switzerland

Uni Graz

Introduction

Recently, a HPTLC-MS interface became available, which semi-automatically can extract zones of interest from TLC/HPTLC plates and direct them into a LC-MS system for substance identification or structure elucidation [1]. So far, the method has hardly been applied for the analysis of plant extracts [2]. Flavonoids are very abundant in plants and play an important role as antioxidants. Therefore, we now have evaluated the application of the HPTLC-MS interface for the investigation of flavonoid containing herbal drugs.



Summary

Rutoside, hyperoside, vitexin, quercetin and rosmarinic acid as pure substances were used to optimize extraction, detection and identification by HPTLC-MS. It was possible to identify hyperoside, vitexin, quercetin and rosmarinic acid also in an extract of *Thymus vulgaris*. The CAMAG HPTLC-MS interface proved to be a quick and powerful tool for the on-line identification of flavonoids in TLC/HPTLC separations. It can complement the classical TLC detection tools and can help to identify individual zones of a TLC/HPTLC separation.

Principle

The instrument is used to locate compounds from a TLC/HPTLC plate and transfer them into a mass spectrometer for identification or structure elucidation. CAMAG TLC-MS interface can be connected to any brand of LC-coupled mass spectrometer.

- Semi-automatic instrument including automatic piston movement for precise seal of the TLC/HPTLC zone on both glass plates and aluminium foil
- Extraction directly from the plate using a suitable solvent delivered by a HPLC pump
- On-line transfer into the mass spectrometer
- Automatic cleaning of the piston between the extractions.




Results and Discussion

Rutoside, hyperoside, vitexin, quercetin and rosmarinic acid could be extracted from TLC/HPTLC silica plates by the CAMAG TLC-MS interface using acetonitrile as solvent delivered by a HPLC pump at 100 µl/min. After application of 2 µg of each compound, significant mass spectra were obtained by electro spray ionization (ESI) in the


negative ion mode. Mass spectrometric signals were obtained within a minute per substance zone. The spectra were suitable for identification of the compounds. It was possible to identify hyperoside, vitexin, quercetin and rosmarinic acid also in an methanolic extract of *Thymus vulgaris*. The applied amount corresponded to 1 mg plant equivalent and was comparable to regular amounts in HPTLC

separations. The CAMAG HPTLC-MS interface proved to be a quick and powerful tool for the on-line identification of flavonoids in TLC/HPTLC separations. It can complement the classical TLC detection tools. The TLC-MS interface can also be used for extraction into vials for NMR, (ATR-)FTIR, static nanospray, direct inlet EI-MS, or MALDI. Scraping off the plate is no more necessary.



Experimental


Extract of *Thymus vulgaris* (500 mg extracted with 5 ml methanol in a 4000 µl vial for 10 min) and pure flavonoids have been applied as bands onto HPTLC plates (Merck) using an automatic TLC sampler ATS 4 (CAMAG). Chromatography was performed in an automatic developing chamber ADC 2 with humidity control (CAMAG). Humidity control: saturated MgCl₂ solution (saturation: 20 ml); developing distance from application position/ lower edge: 62/70; developing solvent: ethyl acetate, formic acid, water (15/1/1); developing time: 18 min; plate drying: 5 min; derivatization reagent: NP-RoseB/PEG. Separated zones were eluted from the plate with the CAMAG TLC-MS interface using acetonitrile as solvent delivered by an HPLC pump (at 100 µl/min). The interface was implemented to a Finnigan LCQ Deca XP Plus ion trap mass spectrometer equipped with an electro spray ionization (ESI) source in negative ion mode.



References


1. Luthmann, H. et al. (2007) *Rap. Commun. Mass Spectrom.* 21: 3772-3776.
2. Reich E, Wöhrner V. (2009) *Chrom. Med.* 75(7): 711-718.

Identification of alkaloids in *Sophora flavescens*




BioLumineX - BioLuminizer

detection: bioluminescence *Vibrio fischeri* 06699




A B C D E F G H I np np np

detection: derivat. anisaldehyde / 366 nm




A - I : chamomile oil

np : 4-nitrophenol



Bio-assays: DPPH

Anti-oxidative properties of essential oils




Silica gel 60 F₂₅₄

Toluene : ethyl acetate
95:5


1: carvacrol methylester; 2: thymol; 3: carvacrol; 4: thyme oil; 5:sage oil; 6:neroli oil; 7: lemon, oil; 8: peppermint oil;
9: rosemary oil; 10: chamomile oil; 11:sweet orange oil; 12: manuka oil; 13: tea tree oil; 14: pine oil; 15: niaouli oil;

Bio-assays: DPPH

Anti-oxidative properties of flavonoids



1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16



1: chlorogenic acid; 2: quercitrin; 3: rutin; 4: ginkgo leaf extract; 5: St. John's wort; 6: great mullein; 7: ribwort plantain; 8: rosemary; 9: majoram; 10: basil; 11: thyme; 12: chamomile; 13: peppermint; 14: arnica; 15: birch; 16: hibiscus

Silica gel 60 F₂₅₄

Ethyl acetate, acetic acid, formic acid, water
100:11:11: 27

55



Thank you!

eike.reich@camag.com

www.camag-laboratory.com

Services

Methods

Publications

Virtual Classroom

About us

CAMAG home

What's new?

New Application Notes

- Content uniformity test of Coenzyme Q10 >>
- HPTLC assay of metronidazole >>
- Determination of apolar lipids from human skin by HPTLC >>
- Determination of ceramides from human skin by HPTLC >>

New Publications

- Quality control of multicomponent herbal drugs: example from TCM >>

HPTLC of dyes



New CAMAG TLC-MS Interface:

For direct extraction of compounds from TLC/HPTLC layers into MS.

[More details >>](#)

CAMAG Laboratory

a center for applied HPTLC

- Broad range of services
- Dedicated to research
- Source of information

Services

Method development and validation

- Download new CAMAG Lab flyer (pdf)
- Feasibility studies
- Contract analyses
- Training

Publications

- Validated HPTLC method for the determination of illegal dyes in spices
- Separation of phospholipids by HPTLC
- Validation of HPTLC identification methods for botanicals

[Download SOP for HPTLC \(pdf\)](#)

[Home](#)

[Contact](#)

Last update: 05.07.2011