

***Planar chromatography meets direct  
ambient mass spectrometry:  
current trends***

**Elizabeth Crawford<sup>1,2</sup>, Brian Musselman<sup>1</sup>**

<sup>1</sup> IonSense, Inc. Saugus, MA, USA

<sup>2</sup> Institute of Chemical Technology, Prague, Czech Republic

# Major Hyphenations of Open Air Ambient MS & (HP)TLC Separations

## Thermal Desorption

- DART – **D**irect **A**nalysis in **R**eal **T**ime
- FAPA – **F**lowing **A**tmospheric **P**ressure **A**fterglow

## Laser Ablation

- LA-DART – **L**aser **A**blation DART
- LAESI – **L**aser **A**blation **E**lectrospray **I**onization

## Liquid Spray

- DESI – **D**esorption **E**lectrospray **I**onization
- EASI – **E**asy **A**mbient **S**onic-spray **I**onization
- LESA – **L**iquid **E**xtraction **S**urface **A**nalysis

- TLC-CAMAG ESI Interface – main ESI based method

# ■ Focusing on *Open Air Thermal Ionization* Coupled with Mass Spectrometry

***Major benefits*** compared to liquid based methods:

- **No use** of solvents during ionization
- **Only semi-destructive** - all material not completely ablated → can run repeat analyses or other methods
- **Complementary MS data** compared with HPTLC ESI-MS – no salt adduct formation, potential thermal separation, both RP & NP plates
  - → **Heating Ramp** (RT – 550° C)
- **Silica gel not disturbed** from plate surface → interaction only with heated gas → **no source contamination**


# DART TLC Evolutions...




**DART-100 Yrs: 2006-09**




**DART-100 & DART-SVP  
Yrs: 2006-present**




**DART-SVP at Angle  
Yrs: 2009-present**






**Laser Ablation DART  
2010 1<sup>st</sup> Publication  
2012 1<sup>st</sup> TLC DART**



# Direct TLC Analysis: Vertical On-Edge Method



**Sample preparation & separation via HPTLC (Methanolic Extract) (Silica gel 60 F<sub>254</sub> at thickness of 200 µm)**




**Helium at 350° C  
0,5 mm/s speed**



**First coupling of DART with TLC:** Morlock, G.; Schwack, W. Determination of isopropylthioxanthone (ITX) in milk, yoghurt and fat by HPTLC-FLD, HPTLC-ESI/MS and HPTLC-DART/MS. *Anal. Bioanal. Chem.* **2006**, 385(3): 586-595.

# Direct analysis vertical edge of HPTLC plate


Benzedrine Tablet: MeOH extract separated by HPTLC (silica gel 60 F<sub>254</sub> at thickness of 200 μm)




# Direct TLC Analysis: *Horizontal Method*




*The height of TLC plate is lower than the He gas flow from ceramic cap !!  
He gas flow and TLC plate must be arranged in a straight line*



Optim. by  
J.P. Yang et  
al. at Kyung  
Hee  
University  
Seoul,  
South Korea




# Normal Phase (NP) Silica Plates (cut narrow)




Data provided by H.J. Kim and Dr. Y.P. Jang (2011) ; College of Pharmacy Kyung Hee University Seoul, South




# Direct TLC Analysis: 45° Angle



A small dab of glycerol was added to the surface of the TLC plate to enhance the heating of the target area on the TLC plate.




The TLC plate was mechanically moved under the DART-SVP source at 0.3 mm/s speed.



Chernetsova, E.S.; Revelsky, A.I.; Morlock, G.E. Some new features of Direct Analysis in Real Time mass spectrometry utilizing the desorption at an angle option. *Rapid Comm. Mass Spec.*, 2011, 25(16): 2275-2282.




# Direct TLC Plate Analysis Workflow

Developed TLC plates are ready for direct mass spectrometric analysis



TLC plate is placed into the holder, which is mounted onto a motorized linear rail

A: 350° C heater setting (Compound 1)  
B: 450° C heater setting (Compound 2)  
C: 450° C plus glycerol on spot (Cmpds 2&3)





# TLC Plate: Pharma Compound 3

Replicate of 10 Repeat Sampling of the Same TLC Spot:  
Glycerol reapplied to single TLC spot before reintroducing the spot to the DART source

450° C

RT: 5.07 - 11.23



Replicates n = 3:  
NO Glycerol applied


# SS1: Well Separated - Lower & Upper Spots

SS1 Plate\_Lower\_Upper Spots\_01  
He; 350C without/with Glycerol, Full MS

8/3/2012 3:14:24 PM


SS1 Plate

RT: 0.00 - 3.47



NL: 7.15E6  
m/z: 233.29-237.68  
MS SS1  
Plate\_Lower\_Upper Spots\_01

NL: 2.42E7  
m/z: 334.27-339.82  
MS SS1  
Plate\_Lower\_Upper Spots\_01




NL: 5.35E5  
SS1 Plate\_Lower\_Upper Spots\_01#510-567 RT: 1.35-1.49 AV: 58 T: ITMS + p NSI Full ms [80.00-1000.00]

NL: 1.05E6  
SS1 Plate\_Lower\_Upper Spots\_01#886-994 RT: 2.31-2.56 AV: 109 T: ITMS + p NSI Full ms [80.00-1000.00]

# For the Next Generations: Laser Ablation TLC Plasma (DART) Ionization

Abstract:

“The PAMLDI-MS system was successfully applied in the detection of low molecular weight compounds from different kinds of samples separated on a normal-phase silica gel, such as dye mixtures, drug standards, and tea extract, with the detection level of 5 ng/mm<sup>2</sup>.”



Zhang, J.; Zhou, Z.; Yang, J.; Zhang, W.; Bai, Y.; Liu, H. Thin Layer Chromatography/Plasma Assisted Multiwavelength Laser Desorption Ionization Mass Spectrometry for Facile Separation and Selective Identification of Low Molecular Weight Compounds. *Anal. Chem.*, **2012**, 84(3): 1496–1503.

# Conclusions

- Both **RP & NP plates** (glass backed) can be directly analyzed without source contamination
- Plates can be cut **thin** (5 mm wide) and run with **DART** source **in-line** with MS
- To analyze intact **large format plates**, set the **DART-SVP** source at an angle (**45° angle**) to the plate
- Resolution of DART gas beam **3-4 mm**, with **laser** target area becomes **nm spot size**
- Coupling with **high resolution accurate mass** (HRAM) MS is major trend with coupling **TLC** and **ambient MS**.
- Direct **quantitation** of HPTLC plates with DART-MS...Tim Häbe & Prof. G. Morlock
  - **O15** today at **14.45**

# DART (HP)TLC Publications

2006

Morlock, G.; Schwack, W. Determination of isopropylthioxanthone (ITX) in milk, yoghurt and fat by HPTLC-FLD, HPTLC-ESI/MS and HPTLC-DART/MS. *Anal. Bioanal. Chem.*, **2006**, 385(3): 586-595.

2007

Morlock, G.; Ueda, Y. J. New coupling of planar chromatography with direct analysis in real time mass spectrometry. *Chromatogr. A*, **2007**, 1143: 243-251.

2008

Dytkiewitz, E. and Morlock, G.E. Analytical strategy for rapid identification and quantification of lubricant additives in mineral oil by high-performance thin-layer chromatography with UV absorption and fluorescence detection combined with mass spectrometry and infrared spectroscopy. *Journal of AOAC International*, **2008**, 91(5): 1237-1243.

Smith, N.J.; Domin, M.A.; Scott, L.T. HRMS Directly From TLC Slides. A Powerful Tool for Rapid Analysis of Organic Mixtures. *Org. Lett.*, **2008**, 10(16): 3493-3496.

2009

Jee, E. H.; Jeong, C. W.; Jeong, S. D.; Kim, H. J.; Jang, Y. P., Detection of characterising compounds on TLC by DART-MS. *Planta Med*, **2009**, 75(9): 38.

2010

Kim, H. J.; Jee, E. H.; Ahn, K. S.; Choi, H. S.; Jang, Y. P. Identification of marker compounds in herbal drugs on TLC with DART-MS. *Archives of Pharmacal Research*, **2010**, 33(9): 1355-1359.

Kim, H. J.; Oh, M. S.; Hong, J.; Jang, Y. P., Quantitative analysis of major dibenzocyclooctane lignans in schisandrae fructus by online TLC-DART-MS. *Phytochemical Analysis*, **2010**, 22(3): 258-262.

2011

Cheng, S.-C.; Huang, M.-Z.; Shiea, J. Thin layer chromatography/mass spectrometry. *Journal of Chromatography A*, **2011**, 1218(19): p. 2700-2711.

Chernetsova, E.S.; Revelsky, A.I.; Morlock, G.E. Some new features of Direct Analysis in Real Time mass spectrometry utilizing the desorption at an angle option. *Rapid Comm. Mass Spec.*, **2011**, 25(16): 2275-2282.

**Howlett, S.E.; Steiner, R.R., Validation of Thin Layer Chromatography with AccuTOF-DART™ Detection for Forensic Drug Analysis. *Journal of Forensic Sciences*, 2011, 56(5): 1261-1267.**

Kim, H.J.; Oh, M.S.; Hong, J.; Jang, Y.P. Quantitative analysis of major dibenzocyclooctane lignans in schisandrae fructus by online TLC-DART-MS. *Phytochemical Analysis*, **2011**, 22(3): 258-262.

**Wood, J.L.; Steiner, R.R. Purification of pharmaceutical preparations using thin-layer chromatography to obtain mass spectra with Direct Analysis in Real Time and accurate mass spectrometry. *Drug Testing and Analysis*, 2011, 3(6): 345-351.**

2012

Morlock, G.E.; Chernetsova, E.S. Coupling of planar chromatography with direct analysis in real time mass spectrometry (DART-MS). *Cent Eur J. Chem.* **2012**, 10(3): 703-710.

Zhang, J.; Zhou, Z.; Yang, J.; Zhang, W.; Bai, Y.; Liu, H. Thin Layer Chromatography/Plasma Assisted Multiwavelength Laser Desorption Ionization Mass Spectrometry for Facile Separation and Selective Identification of Low Molecular Weight Compounds. *Anal. Chem.*, **2012**, 84(3): 1496-1503.

2013

Djelal, H.; Cornée, C.; Tartivel, R.; Lavastre, O.; Abdeltif, A. The use of HPTLC and Direct Analysis in Real Time-Of-Flight Mass Spectrometry (DART-TOF-MS) for rapid analysis of degradation by oxidation and sonication of an azo dye. *Arabian Journal of Chemistry*, **2013**, (0).

# Questions?

Are your spots 'priceless'?



Email: [crawford@ionsense.com](mailto:crawford@ionsense.com)