



University of Hohenheim
Institute of Food Chemistry

Planar SPE coupled to flow injection TOFMS analysis – a rapid pesticide screening tool


Wolfgang Schwack, Claudia Oellig
University of Hohenheim



Planar solid phase extraction

- a new cleanup concept in multi-residue analysis of pesticides


- Application of QuEChERS extracts (50 μ L) onto TLC amino plates
- Twofold development






Planar solid phase extraction

- Application of QuEChERS extracts (50 μ L) onto TLC amino plates
- Twofold development





TLC-LC/MS



Extraction of target zones (acetonitrile/ammonium formate)

- on-line coupling to LC/MS column (Chromolith RP-18)
- elution into autosampler vials (LC/MS or GC/MS)

→ Very clean, matrix-free extracts, no matrix effects!



pSPE- μ L-FIA-TOFMS

University of Hohenheim
Institute of Food Chemistry

pSPE





TLC-MS



μ L-FIA



nano-spray TOFMS





pSPE- μ L-FIA-TOFMS

University of Hohenheim
Institute of Food Chemistry

pSPE





TLC-MS



μ L-FIA




nano-spray TOFMS



- No chromatography, all pesticides in a single FIA peak.
- Whole sample information in a single mass spectrum.





Example cucumber extracts






Example cucumber extracts







Zoom





Target / Non-target screening



EXCEL Sheet 1: data base table (about 250 pesticides)

EXCEL Sheet 2: results table (TOFMS mass list)

ACCESS Query: comparison of tables

- positive correlation = target screening

- negative correlation = non-target screening



Real samples - pSPE-FIA-TOFMS ↔ dSPE-LC-MS/MS ¹⁾

| Sample | Pesticide | LC-MS/MS [mg/kg] | FIA-TOFMS [S/N] |
|----------------------|--------------------|------------------|-----------------|
| Banana | Thiabendazole | 0.76 | 2004 |
| | Imazalil | 0.73 | 1090 |
| Blackberry | Cyprodinil | 0.86 | 2347 |
| | Lambda-cyhalothrin | 0.01 | n.d. |
| | Thiacloprid | 0.01 | n.d. |
| Currant | Cyprodinil | 0.47 | 1373 |
| | Boscalid | 0.37 | 124 |
| | Pyraclostrobin | 0.14 | 121 |
| | Trifloxystrobin | 0.06 | n.d. |
| Savoy cabbage | Azoxystrobin | 1.3 | 1770 |
| | Difenoconazole | 0.39 | 317 |
| | Pymetrozine | 0.20 | n.d. |
| | Cyfluthrin | 0.09 | n.d. |
| | Indoxacarb | 0.03 | n.d. |

1) Institute Kirchhoff, Berlin



Conclusions

- Even our low-resolution TOFMS provided exact results.
- The novel pSPE-FIA-TOFMS concept fully meets the SANCO document criteria for a qualitative screening method.
- FIA-TOFMS clearly showed the high efficiency of planar SPE.
- All residues are detected in a single peak/mass spectrum.
- A simple and rapid target screening is allowed.
- A true non-target screening is also simply enabled.
- A state-of-the art HRMS undoubtedly will improve the hit rates.


👉 [Poster P-37 \(experimental details\)](#)



University of Hohenheim
Institute of Food Chemistry





Many thanks!

Claudia Oellig





Matrix effects \leftrightarrow recoveries (n=5)



post-extraction spiking